# **Energy sensors**



# 1.2 Energy Sensors

### Introduction

Pyroelectric sensors are for measuring repetitive pulse energies and average powers at pulse rates up to 25000 pulses per second and pulse widths up to 20ms. Note that single shot energy with pulse rates less than one pulse every 5s or so can be measured with thermal sensors described in the power sensor section.

### **Pyroelectric Sensors**

Pyroelectric type sensors are useful for measuring the energy of repetitively pulsed lasers at up to 25,000Hz and are sensitive to low energies.

They are less durable than thermal types and therefore should not be used whenever it is not necessary to measure the energy of each pulse and average power measurement is sufficient.


Pyroelectric sensors use a pyroelectric crystal that generates an electric charge proportional to the heat absorbed. Since the two surfaces of the crystal are metalized, the total charge generated is collected and therefore the response is not dependent on beam size or position. This charge then charges a capacitor in parallel with the crystal and the voltage difference thus generated is proportional to the pulse energy. After the energy is read by the electronic circuit, the charge on the crystal is discharged to be ready for the next pulse. The response time of the pyroelectric sensor depends on the time it takes for the heat to enter the crystal and heat it up. For metallic type pyro detectors, this time is tens of  $\mu$ s and thus the metallic type can run at a high repetition rate. For the BF and BB type, the response time is hundreds of  $\mu$ s with a correspondingly lower repetition rate.

Ophir pyroelectric detectors have unique and proprietary circuitry that allow them to measure long pulses as well as short pulses and work at a high duty cycle, i.e. where the pulse width is as much as 30% of the total cycle time.

Ophir came out with the compact C line of pyroelectric sensors that replaced previous models. The electronics and mechanics has been completely upgraded and the current sensors are superior in every way: more compact, wider dynamic range, have higher repetition rates and measure longer pulses. Through constant development, Ophir again brings you the best performance in the market.


Note: Older line of Pyroelectric sensors is not supported by the StarBright and StarLite meters.

All Ophir power and energy sensors come with a mounting stand.




# **Absorption and Damage Graphs for Pyroelectric Sensors**

### Absorption vs. Wavelength




# Damage Threshold vs. Pulse Width




# **Wavelength Range and Repetition Rate for Energy Sensors**

### **Wavelength Range**



### **Repetition Rate Range**



# 1.2.1 Photodiode Energy Sensors

### 10pJ to 15μJ

#### **Features**

- Silicon and Germanium detectors
- Very sensitive down to 10pJ
- Repetition rates to 20kHz
- Wide spectral range

### PD10-C / PD10-IR-C / PD10-pJ-C / PD10-IR-pJ-C




Model	PD10-C		PD10-IR-C	PD10-IR-C			PD10-IR-pJ-C		
Use	Low energies	s	Infrared		Lowest ener	gies	Infrared, I	lowest energies	
Aperture mm	Ø10		Ø5		Ø10		Ø5		
Absorber Type	Si photodiode	ة	Ge photodiode	ة	Si photodiod	e	Ge photodiode		
Spectral Range µm (a)	0.19 - 1.1		0.7 – 1.8		0.2 - 1.1		0.7 - 1.8		
Surface Reflectivity % approx.	50		30		30		30		
Calibration Accuracy +/-% (a)	5		5		5		5		
Energy Scales	20µJ to 20nJ		600nJ to 6nJ		200nJ to 200	pJ	20nJ to 20	0pJ	
Lowest Measurable Energy nJ (b)	1 at 900nm		1 at 1550nm		0.01 at 900nr	n	0.03 at 155	50nm	
Max Pulse Width ms	0.005		0.005		0.005		0.005		
Maximum Pulse Rate pps	20kHz		10kHz		20kHz		10kHz		
Noise on Lowest Range nJ	0.05		0.1		0.001		0.01		
Additional Error with Frequency %	±1% to 20kHz	, (c)	±1.5% to 10kH:	Z	±1% to 20kH	Z <sup>(d)</sup>	±1.5% to 1	0kHz	
Linearity with Energy for > 10% of full scale (b)	±1.5%		±1.5%		±1.5%		±1.5%		
Damage Threshold J/cm <sup>2</sup>	0.1		0.1		0.1		0.1		
Maximum Average Power mW	50 at 800nm		6		0.5		0.1		
Maximum Average Power Density			50		5		5		
W/cm <sup>2</sup>	30		30		5		5		
Maximum Energy vs. Wavelength	Wavelength	Max Energy	Wavelength	Max Energy	Wavelength	Max Energy	Wavelengtl	h Max Energy	
	<300nm	15μJ	800 - 900nm	600nJ	<300nm	150nJ	800 - 900nr	m 20nJ	
	350 - 550nm	8µJ	1000 - 1300nm	200nJ	350 - 550nm	75nJ	1000 - 1300	Onm 8nJ	
	>800nm	5μJ	1300 - 1400nm	170nJ	>800nm	50nJ	1300 - 1400	Onm 7nJ	
			1480 - 1560nm	ı 150nJ			1480 - 156	i0nm 6nJ	
			>1650nm	600nJ			>1650nm	20nJ	
Fiber Adapters Available (see page 102)	ST, FC, SMA, S	C	ST, FC, SMA, SC		ST, FC, SMA, S	SC .	ST, FC, SM	A, SC	
Weight kg	0.25		0.25		0.25		0.25		
Compliance	CE, China Rol-	lS .	CE, China RoHS	5	CE, China Rol	HS	CE, China I	RoHS	
Version	,		• •		•		·		
Part number	7Z02944		7Z02955		7Z02945		7Z02946		
Note: (a) This is basic calibration accuracy. In certain wavelength regions calibration there is additional error as tabulated here.		d ±3% d ±2%		±2% ±2%		ld ±2% ld ±2%	<900nm >1700nm	add ±2% add ±2%	

Note: (b) With the "user threshold" setting set to minimum. For other settings, the spec is for >10% of full scale or greater than twice the "user threshold", whichever is greater. The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PD-C series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments. For further information, see the FAQs on our Website.


Note: (c) Additional Error with Frequency of  $\pm 1\%$  on ly for energies up to  $2\mu J$ . For higher energies  $\pm 1\%$  up to 10kHz, -4% at 20kHz.

Note: (d) Additional Error with Frequency of  $\pm 1\%$  only for energies up to 20nJ. For higher energies  $\pm 2\%$  up to 10kHz, -5% at 20kHz.

### PD10-C / PD10-pJ-C



### PD10-IR-C / PD10-IR-pJ-C





### 0.1μJ to 1mJ

#### **Features**

- Ø8mm aperture
- Repetition rates up to 20,000Hz
- High sensitivity sensors
- Pulse widths up to 20µs





Model	PE9-C			PE9-ES-C		
Use	Very Sensitiv	e		Most Sensitive		
Aperture mm	Ø8			Ø8		
Absorber Type Spectral Range µm <sup>(a)</sup>	metallic			metallic		
Spectral Range µm (a)	0.15 - 12			0.15 - 12		
Surface Reflectivity % approx.	50			50		
Calibration Accuracy +/-% (a)	3			3		
Max Pulse Width Setting (c)	1µs	2µs	20µs	1µs	2µs	20µs
Energy Scales	1mJ to 2µJ	1mJ to 2µJ	1mJ to 20µJ	200µJ to 200nJ	200µJ to 200nJ	200µJ to 2µJ
Lowest Measurable Energy µJ (b)	0.5	< 0.2	0.5	0.1	<0.1	0.1
Max Pulse Width µs	1	2	20	1	2	20
Maximum Pulse Rate pps	25kHz	15kHz	10kHz	20kHz	15kHz	10kHz
Noise on Lowest Range µJ	0.04	0.05	0.1	0.01	0.01	0.02
Additional Error with Frequency %	±1% to 15kHz, ±6% to 25kHz	±1% to 15kHz	±1% to 10kHz	±1.5% to 20kHz	±1.5% to 15kHz	±1.5% to 10kHz
Damage Threshold J/cm <sup>2</sup>						
<100ns	0.1			0.1		
1µs	0.2			0.2		
300µs	3			3		
Linearity with Energy (b)	±1%			±1.5%		
Maximum Average Power W	2			2		
Maximum Average Power Density W/cm <sup>2</sup>	30			30		
Fiber Adapters Available (see page 102)	ST, FC, SMA, SC	_		ST, FC, SMA, SC		
Weight kg	0.25			0.25		
Compliance	CE, China RoH	S		CE, China RoHS		
Version	,			,		
Part Number	7Z02933			7Z02949		
Note: (a) Calibrated curve is checked and adjusted at the following wavelengths (µm)	0.193, 0.355, 1.064,	1.48-1.6		0.355, 1.064, 1.48-1.6		
For other wavelengths in the curve there is additional calibration error as stated.	240-800nm add ±44	%, 2-3μm add ±8%, 1	0.6μm add ±15%.	240-800nm add ±4%, 2 <240nm not calibrated	2-3μm add ±8%, 10.6μm a	add ±15%.


Note: (b) For >7% (>10% for PE9-ES-C) of full scale, with the "user threshold" setting set to minimum. For other settings, the spec is for >7%/>10% of full scale or greater than twice the "user threshold", whichever is greater.

The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PE-C series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments.

For further information, see the FAQs on our Website.

Note: (c) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter, only 2 out of 3 pulse widths settings are available; the 1µs (displayed as "10µs") and the 2µs (displayed as "20µs").

#### PE9-C / PE9-ES-C





### 1μJ to 10mJ

#### **Features**

- Ø12mm apertures
- Repetition rates up to 25,000Hz
- High sensitivity sensors
- Pulse widths up to 5ms

#### PE10-C / PE10BF-C



Model	PE10-C		PE10BF-C	
Use	Sensitive		High damage threshold	
Aperture mm	Ø12		Ø12	
Absorber Type	metallic		BF	
Spectral Range µm (a)	0.15 - 12		0.15 - 3, 10.6 <sup>(d)</sup>	
Surface Reflectivity % approx.	50		20	
Calibration Accuracy +/-% (a)	4		3	
Max Pulse Width Setting (e)	1µs	30µs	1ms	5ms
Energy Scales	10mJ to 2µJ	10mJ to 20μJ	10mJ to 20µJ	10mJ to 200µJ
Lowest Measurable Energy µJ (c)	1	1	7	20
Max Pulse Width μs	1	30	1000	5000
Maximum Pulse Rate pps	25kHz	5kHz	250Hz	50Hz
Noise on Lowest Range µJ	0.1	0.15	1	5
Additional Error with Frequency %	±2% to 15kHz, ±3% to 25kHz	±1% to 5kHz	±1%	±1%
Damage Threshold J/cm <sup>2</sup>				
<100ns	0.1		0.8 <sup>(b)</sup>	
1µs	0.2		1 <sup>(b)</sup>	
300µs	3		4 <sup>(b)</sup>	
Linearity with Energy (c)	±1.5%		±2%	
Maximum Average Power W	2		3	
Maximum Average Power Density W/cm <sup>2</sup>	50		50	
Fiber Adapters Available (see page 102)	ST, FC, SMA, SC		ST, FC, SMA, SC	
Weight kg	0.25		0.25	
Compliance	CE, China RoHS		CE, China RoHS	
Version				
Part Number	7Z02932		7Z02938	
Note: (a) Calibrated curve is checked and adjusted at the following wavelengths (μm)	1.064, 0.355		0.193, 0.248, 0.355, 0.532, 1.064	
For other wavelengths in the curve there is additional calibration error as stated.	240 - 800nm add ±4%, 2-3µm a <240nm not calibrated	aα ±8%, 10.6μm add ±15%.	0.2-3μm ±2%, 10.6μm ±5%	


Note: (b) For wavelenghts below 600nm, derate damage threshold to 60% of given values. Below 300nm, derate to 40% of given values.

Note: (c) For >7% of full scale, with the "user threshold" setting set to minimum. For other settings, the spec is for >7% of full scale or greater than twice the "user threshold", whichever is greater. The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PE-C series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments. For further information, see the FAQs on our Website.

Note: (d) The absorption at 675nm is approximately the same as at 10.6µm. Therefore, to measure a CO2 laser, set to the 675nm setting. The additional error for measuring 10.6µm is ±5%.

Note: (e) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter, for the PE10-C model the 1µs pulse width setting is displayed as "10µs".

#### PE10-C / PE10BF-C



### 8μJ to 10J

#### **Features**

- Ø24mm apertures
- Metallic coating for high rep rates
- BF coating for highest damage threshold
- Rep rates up to 10kHz
- Measure lasers with pulse widths up to 20ms





PE25BF-C



Model	PE25-C		PE25BF-C								
Use	High rep	rate				High damage threshold					
Aperture mm	Ø24					Ø24					
Absorber Type	metallic					BF					
Spectral Range µm (a)	0.15 - 3					0.15 - 3,	10.6 (e)				
Surface Reflectivity % approx.	50					20					
Calibration Accuracy +/-% (a)	3					3					
Max Pulse Width Setting (d)	2µs	30µs	500µs	1ms	5ms	1ms	2ms	5ms	10ms	20ms	
Energy Scales	10J to	10J to	10J to	10J to	10J to	10J to	10J to	10J to	10J to	10J to	
- 3,	200uJ	200µJ	2mJ	2mJ	2mJ	2mJ	2mJ	20mJ	20m J	20mJ	
Lowest Measurable Energy µJ (c)	8	10	60	80	100	60	100	400	400	400	
Max Pulse Width ms	0.002	0.03	0.5	1	5	1	2	5	10	20	
Maximum Pulse Rate pps	10kHz	5kHz	900Hz	450Hz	100Hz	250Hz	100Hz	50Hz	40Hz	20Hz	
Noise on Lowest Range µJ	0.5	1	6	10	20	10	20	40	40	50	
Additional Error with Frequency %	+2% to	+1.5%	±2% to		±1.5% to		+1%	+1%	+1%	±2%	
Additional Error with rrequerity 70	5kHz ±4% to 10kHz		750Hz	400Hz	80Hz	1170	11/0	1170	±170	1270	
Linearity with Energy for >7% of full scale (c)	±1.5%					±2%					
Damage Threshold J/cm <sup>2</sup> (b)											
<100ns	0.1					0.8					
1µs	0.2					1					
300µs	2					4					
2ms	6					10					
Maximum Average Power W	15, 25 with	n optional	heat sink			15, 25 wi	ith optiona	al heat sink	(		
Maximum Average Power Density W/cm <sup>2</sup>	20	'				20					
Uniformity over surface	±2% over	r central :	50% of ap	erture		±2% ov	er central	150% of a	perture		
Fiber Adapters Available (see page 102)	ST, FC, SN					ST, FC, S			•		
Weight kg	0.25	,				0.25	,				
Compliance	CE, China	RoHS				CE, Chir	na RoHS				
Version											
Part Number	7Z02937	,				7Z0293	35				
Note: (a) Calibration curve is verified and adjusted at specified wavelengths.	Specified wavelengths: 248-266nm, 355nm, 1064nm and Specifie						Specified wavelengths: 193nm, 248 266nm, 355nm, 532nm and 1064nm.				
At other wavelengths, there may be an additional error up to the value given.								6.			
Note: (b)	For wavelengths below 600						0nm, derate damage threshold				
						to 60% of given valu		s. Below 300	nm, derate t	o 40% of	

Note: (c) With the "user threshold" setting set to minimum. For other settings, the spec is for >7% of full scale or greater than twice the "user threshold", whichever is greater. The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PE-C series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments. For further information, see the FAQs on our Website.

Note: (d) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter, only 2 out of 5 pulse widths settings are available; for the PE25-C model the 2µs (displayed as "10µs") and 1ms settings, and for the PE25BF-C model the 1ms and 10ms settings.

Note: (e) If the sensor is set to the 1064nm wavelength, then when measuring 10.6µm pulses, the reading will be approximately 1.19X the correct reading. If you use the attenuate function and set the attenuation to read 0.84, then you will have the correct reading at 10.6µm. The additional error at 10.6µm is +/-5%.



<sup>\*</sup> For drawings please see page 99

### 10μJ to 10J

#### **Features**

- Ø46mm apertures
- Metallic coating for high rep rates
- BF coating for highest damage threshold
- Rep rates up to 10kHz
- Measure lasers with pulse widths up to 20ms



PE50-C



PE50BF-C





Model	PE50-C					PE50BF-0	c					
Use	High rep	rate				High damage threshold						
Aperture mm	Ø46					Ø46						
Absorber Type	metallic					BF						
Spectral Range µm (a)	0.15 - 3					0.15 - 3, 1	0.6 <sup>(e)</sup>					
Surface Reflectivity % approx.	50					20						
Calibration Accuracy +/-% (a)	3					3						
Max Pulse Width Setting (d)	2µs	30µs	500µs	1ms	5ms	1ms	2ms	5ms	10ms	20ms		
Energy Scales	10J to 200µJ	10J to 200µJ	10J to 2mJ	10J to 2mJ	10J to 2mJ	10J to 2mJ	10J to 2mJ	10J to 20mJ	10J to 20mJ	10J to 20mJ		
Lowest Measurable Energy µJ (c)	10	10	60	80	100	120	300	600	600	600		
Max Pulse Width ms	0.002	0.03	0.5	1	5	1	2	5	10	20		
Maximum Pulse Rate pps	10kHz	5kHz	900Hz	450Hz	100Hz	250Hz	100Hz	50Hz	40Hz	20Hz		
Noise on Lowest Range µJ	0.5	1	6	10	20	30	60	100	100	100		
Additional Error with Frequency %	±2% to 2kHz ±4.5% to 5kHz	±2%	±2% to 750Hz	±2% to 400Hz	±1% to 80Hz	±1%	±1%	±1%	±1%	±2%		
Linearity with Energy for >7% of full scale (c)	±1.5%					±2%						
Damage Threshold J/cm <sup>2 (b)</sup>	_11.570											
<100ns	0.1					0.8						
1µs	0.2					1						
300µs	2					4						
2ms	6					10						
Maximum Average Power W	15. 25 with	optional he	at sink			15. 25 with	optional h	eat sink				
Maximum Average Power Density W/cm <sup>2</sup>	20					20						
Uniformity over surface	+2% over	central 509	% of apertu	ıre		+2% over	central 50	% of apertu	ıre			
Fiber Adapters Available (see page 102)	ST, FC, SM					ST, FC, SM		,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,				
Weight kg	0.25	,				0.25	,					
Compliance	CE, China	RoHS				CE, China	RoHS					
Version	,											
Part Number	7Z02936					7Z02934						
Note: (a) Calibration curve is verified and adjusted at specified wavelengths. At other wavelengths, there may be an additional error up to the value given.	Specified wavelengths: 248-266nm, 355nm and 1064nm.  Max additional error at 2940nm ±3%.  Max additional error at other wavelengths: ±2%.  <240nm not calibrated					Specified wavelengths: 193nm, 248-266nm, 355nm, 532nm and 1064nm.  Max additional error at 2940nm ±3%.  Max additional error at other wavelengths: ±2%.						
Note: (b)	For wavelengths below 600nm, derate damage thresh of given values. Below 300nm, derate to 40% of given v											

Note: (c) With the "user threshold" setting set to minimum. For other settings, the spec is for >7% of full scale or greater than twice the "user threshold", whichever is greater. The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PEC series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments. For further information, see the FAQs on our Website.

Note: (d) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter, only 2 out of 5 pulse widths settings are available; for the PE50-C model the 2µs (displayed as "10µs") and 1ms settings, and for the PE50BF-C model the 1ms and 10ms settings.

Note: (e) If the sensor is set to the 1064nm wavelength, then when measuring 10.6µm pulses, the reading will be approximately 1.19X the correct reading. If you use the attenuate function and set the attenuation to read 0.84, then you will have the correct reading at 10.6µm. The additional error at 10.6µm is +/-5%.



<sup>\*</sup> For drawings please see page 99

### 20μJ to 10J

#### **Features**

- Sensors with diffuser for high energies and high energy densities
- Metallic coating for high rep rates
- BF coating for highest damage threshold
- Wide spectral range. Measure YAG and harmonics and many more.
- Rep rates up to 10kHz
- Measure lasers with pulse widths up to 20ms



Model	PE50-D	IF-C				PE25BF-DIF-C						
Use	High re	p rate. Co	mplete ca	libration	curve	Comple thresho	te calibra Id	tion curv	e. High da	mage		
Aperture mm	Ø35					Ø20						
Absorber Type	Metallic	with diffu	ser			BF with diffuser						
Spectral Range µm (a)	0.19 - 2.2	2, 2.94				0.24 - 2.2	2					
Surface Reflectivity % approx.	25					25						
Calibration Accuracy +/-% (a)	3					3						
Max Pulse Width Setting (d)	2µs	30µs	500µs	1ms	5ms	1ms	2ms	5ms	10ms	20ms		
Energy Scales	10J to 200µJ	10J to 200µJ	10J to 2mJ	10J to 2mJ	10J to 20mJ	10J to 2mJ	10J to 2mJ	10J to 20mJ	10J to 20mJ	10J to 20mJ		
Lowest Measurable Energy µJ (c)	20	20	100	120	200	100	150	200	200	300		
Max Pulse Width ms	0.002	0.03	0.5	1	5	1	2	5	10	20		
Maximum Pulse Rate pps	10kHz	5kHz	900Hz	450Hz	100Hz	250Hz	100Hz	50Hz	40Hz	20Hz		
Noise on Lowest Range µJ	1	2	20	20	40	15	30	40	40	60		
Additional Error with Frequency %	±2% to 2kHz ±4.5% to 5kHz	±2%	±1% to 750Hz	±2% to 400Hz	±1% to 80Hz	±1%	±1%	±1%	±1%	±2%		
Linearity with Energy for >7% of full scale (c)	±1.5%					±2%						
Damage Threshold J/cm <sup>2 (b)</sup>	_11570											
<100ns	1					4						
1µs	2					5						
300µs	20					20						
2ms	40					60						
Maximum Average Power W		ith ontion	al heat sin	k			ith option	al heat sin	k			
Maximum Average Power Density W/cm <sup>2</sup>	100	ппораоп	ai ricat sii i			120	itii optioiii	ai ricat siri	IX.			
Uniformity over surface		er central 2	20mm				er central 1	0mm				
Weight kg	0.25	rei ceritiai 2	20111111			0.25	ci ccittiai i	OTTITT				
Compliance	CE, Chin	a RoHS				CE, Chin	a RoHS					
Version	CL, CIIII	ia non is				CL, CIIII	u 110115					
Part Number	7Z0293	9				7Z0294	1					
Notes: (a) Calibration curve is verified and adjusted at		wavelengths	S:				wavelengths	:				
specified wavelengths.	193nm, 248-266nm, 532nm, 1064nm and 2100nm. 248-266nm, 355nm, 532nm, 1064nm and 2						m and 2100r	nm.				
At other wavelengths, there may be an additional error up to the value given.	Max addit above: ±2	ional error at %.	t 193nm ±4% t other wavel	engths not s		Max additional error at other wavelengths not specified above: $\pm 2\%$ .						
Notes: (b)	193nm reading may need 1min irradiation to stabilize.  For wavelengths > 2.1 µm, derate to 40% of above values.  For beam size <=5mm. For 10mm beam, derate to 40% of above value.  For wavelengths below 600nm, derate to 60% of giver For wavelengths below 240nm, derate to 1J/cm².  For beam size <=4mm. For 8mm beam, derate to 50% above values.							2.				

Notes: (c) With the "user threshold" setting set to minimum. For other settings, the spec is for >7% of full scale or greater than twice the "user threshold", whichever is greater. The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PE-C series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments. For further information, see the FAQs on our Website.

Notes: (d) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter, only 2 out of 5 pulse widths settings are available; for the PE50-DIF-C model the 2µs (displayed as "30µs") and 1ms settings, and for the PE25BF-DIF-C model the 1ms and 10ms settings.



<sup>\*</sup> For drawings please see page 99

### 100µJ to 40J

#### **Features**

- Sensors with diffuser for high energies and high energy densities
- BF coating for highest damage threshold
- BB coating for spectral flatness
- Wide spectral range. Measure YAG and harmonics and many more.
- Rep rates up to 250Hz
- Measure lasers with pulse widths up to 20ms
- PE50BF-DIFH-C sensor highest damage threshold

#### PE50BF-DIF-C / PE50BF-DIFH-C



#### PE50BB-DIF-C

DIFFUSER IN DIFFUSER OUT





Model	PE50BF-DIF-C / PE50BF-DIFH-C					PE50BB-DIF-C							
Use	Comple thresho	ete calibrat old	ion c	urve	e. Highest	damage	Remov	vable di	ffuser. Spe	pectrally flat			
Diffuser	Fixed						Diffuse	r out		Diffuser in			
Aperture mm	Ø35						Ø46			Ø33			
Absorber Type	BF with	diffuser					BB			BB with diffuser			
Spectral Range µm (a)	0.19 - 2.	2, 2.94					0.19 - 1	20		0.4 – 2.5			
Surface Reflectivity % approx.	25						5			15			
Calibration Accuracy +/-% (a)	3						3			3			
Max Pulse Width Setting (d)	1ms	2ms	5m	S	10ms	20ms	3ms	10ms	20ms	3ms	10ms	20ms	
Energy Scales	10J to	10J to	10J	to	10J to	10J to	10J to	10J to	10J to	40J to	40J to	40J to	
37	2mJ	2mJ	20m	٦J	20mJ	20mJ	2mJ	20mJ	20mJ	8mJ	8mJ	8mJ	
Lowest Measurable Energy mJ (c)	0.2	0.4	0.8		0.8	0.8	0.1	0.1	0.2	0.5	5	5	
Max Pulse Width ms	1	2	5		10	20	3	10	20	3	10	20	
Maximum Pulse Rate pps	250Hz	100Hz	50H	lz	40Hz	20Hz	40Hz	10Hz	5Hz	40Hz	10Hz	5Hz	
Noise on Lowest Range µJ	40	80	200		200	200	15	15	20	40	60	80	
Additional Error with Frequency %	±1%	±1%	±1%		±2%	±2%	±1%	+1%	±1%	±1%	±1%	±1%	
Linearity with Energy for >7% of full scale (c)	±2%					±2%							
Damage Threshold J/cm <sup>2</sup> (b)	PE50BF-DIF-C			PE:	50BF-DIFH-	-C	Diffuse	er out		Diffuse	r in		
<100ns	4			6			0.3			3			
1µs	5			8			0.3			3			
300us	20			30			1			10			
2ms	60			90			2			20			
Maximum Average Power W		ith optiona	ıl heat	t sink	k			with opt	ional heat	eat 30, 50 with optional heat sink			
Maximum Average Power Density W/cm <sup>2</sup>	200						10			500			
Uniformity over surface		er central 20	0mm					er 70% o	f diameter		ver centra	l 20mm	
Weight kg	0.25						0.25						
Compliance	CE, Chin	a RoHS		CF	, China Rol	4S	CF Chi	ina RoHS					
Version	,			,	,								
Part Number	7Z0294	0		7Z	02943		7Z029	47					
Notes: (a) Calibration accuracy at various wavelengths as specified here.	Specified 355nm, 53 2100nm.	wavelengths: 32nm, 1064nn		Spe 193 532	ecified wavele 3nm, 248-266 2nm, 1064nm	nm, 355nm, and 2100nm.	Calibrated at 1064nm			Calibrated at 1064nm, 532nm and 2100nm only. Calibration accuracy at 2100nm, ±5%.			
At other wavelengths, there may be an additional error up to the value given.					litional erro gths is ±29								
Notes: (b)	For wavelengths >2.1µm, derate to 10% of above values. For wavelengths below 600nm, derate to 60% of given values (for DIFH 50% of given values). For wavelengths below 240nm, derate to 1J/cm². For beam size <=5mm. For 10mm beam, derate DIF to 80% and DIFH to 70% of above.												


Notes: (c) With the "user threshold" setting set to minimum. For other settings, the spec is for >7% of full scale or greater than twice the "user threshold", whichever is greater. The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PE-C series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments. For further information, see the FAQs on our Website.

Notes: (d) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter only 2 of the pulse width settings are available. For the PE-BF models the 1ms and 10ms settings and for the PE-BB model the 3ms and 10ms settings. Furthermore, with the diffuser mounted, the sensor may saturate at lower than the maximum energy in some cases. Therefore it is recommended to use these sensors with the newer meters/PC interfaces.




<sup>\*</sup> For drawings please see page 99


### PE25-C / PE25BF-C




#### PE25BF-DIF-C




### PE50BF-DIFH-C




### PE50-C / PE50BF-C



#### PE50BF-DIF-C / PE50-DIF-C



### PE50BB-DIF-C



### 10μJ to 40J

#### **Features**

- Removable diffusers
- PE50-DIF-ER-C mainly for NIR lasers
- PE100BF-DIF-C for very large beams
- Rep rates up to 10kHz
- Measure lasers with pulse widths up to 20ms

#### PE50-DIF-ER-C

#### DIFFUSER IN DIFFUSER OUT





#### PE100BF-DIF-C

DIFFUSER IN DIFFUSER OUT






Model	PE50	-DIF-I	R-C								PE100BF-DIF-C										
Use	Main	ly for	1064r	ım, 2.	lμm a	nd 2.9	4μm				Very large aperture										
Diffuser	Diffu:	ser ou	t			Diffus	ser in				Diffu	ser out	t			Diffu	ser in				
Aperture mm	Ø46					Ø33					Ø96					Ø85					
Absorber Type	Meta	llic				Meta	llic wit	h diffu	iser		BF				BF with diffuser						
Spectral Range µm (a)	0.19 -	- 3				0.4 - 3	3				0.15 - 3				0.4 - 2.5						
Surface Reflectivity % approx.	50					50					20					50					
Calibration Accuracy +/-% (a)	3					4					3					4					
Max Pulse Width Setting (c) Energy Scales			500µs 10J to 2mJ			2μs 30J to 600μJ	30J to				1ms 10J to 2mJ	2ms 10J to 20mJ		10ms 10J to 20mJ		40J to	2ms 40J to 40mJ		40J to		
Lowest Measurable Energy mJ (b)	0.01	0.01		0.08	0.1			0.3	0.4	0.5	0.4	0.7			1.5	2	3	5	5	5	
Max Pulse Width ms	0.002	0.03	0.5	1	5	0.002	0.03	0.5	1	5	1	2	5	10	20	1	2	5	10	20	
Maximum Pulse Rate pps	10kHz	5kHz	800Hz	400Hz		10kHz	5kHz			100Hz		100	50	35	25	200	100	50	35	25	
Noise on Lowest Range µJ	1	1	6	10	20	5	5	30	50	100	80	150	250	200	200	300	500	1000	600	600	
Additional Error with Frequency %	2kHz ±4.5% to 5kHz	4.5% 80Hz ±4.5% 80Hz				to															
Linearity with Energy for > 10% of full scale (b)					±1.	.5%									±	1%					
Damage Threshold J/cm <sup>2</sup>																					
<100ns	0.1					1.5					8.0					3					
1µs	0.2					3					1					3					
300µs	2					20					5					10					
2ms	6					60					10					25					
Maximum Average Power W		5 with	option	nal hea	t sink	40, 60	) with	optior	nal hea	t sink	25					50					
Maximum Average Power Density W/cm <sup>2</sup>	20					500					20					500					
Weight kg	0.3										1.2										
Compliance	CE, C	hina R	oHS								CE, C	hina R	oHS								
Version																					
Part Number	7Z02948					7Z02942															
Notes: (a)	Calibrated at 532nm and 1064nm Calibrated at 1064nm, only 2100nm and 2940nm			Calibrated at 532nm and 1064nm only			Calibrated at 532nm, 1064nm and 1550nm only														


Notes: (b) With the "user threshold" setting set to minimum. For other settings, the spec is for >10% of full scale or greater than twice the "user threshold", whichever is greater. For use with Centauri, StarBright, StarLite, Nova II, Vega, Juno, Juno+ and EA-1. The sensors will operate with older Ophir meters and PC interfaces but do not support the threshold function and may give inaccurate readings with the diffuser in and therefore it is not recommended to use these sensors with older Ophir meters and PC interfaces. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments. For further information, see the FAQs on our Website.

Notes: (c) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter only 2 of the 5 pulse width settings are available. For the PE50-DIF-ER-C, the 30µs and 1ms settings and for the PE100BF-DIF-C, the 1ms and 10ms settings. Furthermore, with the diffuser mounted, the sensor may saturate at lower than the maximum energy in some cases. Therefore it is recommended to use these sensors with the newer meters/PC interfaces.

### PE50-DIF-ER-C



#### PE100BF-DIF-C



### 1mJ to 40J

#### **Features**

- Fan or conduction cooled for high average power capability
- BF coating with diffuser for highest damage threshold
- Wide spectral range. Measure YAG and harmonics and many more
- Rep rates up to 250Hz
- Measure lasers with pulse widths up to 20ms





Model	FPE80BF-DIF-C					PE80BF-DIF-C						
Use	High a	verage p	ower pu	lsed laser	's	Large a	perture	pulsed la	sers			
Diffuser	Fixed					Fixed						
Aperture mm	Ø53											
Absorber Type	BF with diffuser BF with diffuser											
Spectral Range µm (a)	0.19 - 2	2.2, 2.94				0.19 – 2.2, 2.94						
Surface Reflectivity % approx.	25					25						
Calibration Accuracy +/-% (a)	3					3						
Max Pulse Width Setting (d)	1ms	2ms	5ms	10ms	20ms	1ms	2ms	5ms	10ms	20ms		
Energy Scales	40J to	40J to	40J to	40J to	40J to	40J to	40J to	40J to	40J to	40J to		
1 (af)	40mJ	40mJ	40mJ	40mJ	40mJ	40mJ	40mJ	40mJ	40mJ	40mJ		
Lowest Measurable Energy mJ (c, f)	1	1	1	2	2	4	4	4	4	4		
Max Pulse Width ms	1	2	5	10	20	1	2	5	10	20		
Maximum Pulse Rate pps	250Hz	100Hz	50Hz	40Hz	20Hz	250Hz	100Hz	50Hz	40Hz	20Hz		
Noise on Lowest Range μJ	200	300	300	300	300	100	200	200	200	200		
Additional Error with Frequency %	±1.5% ±1.5% ±1.5% ±1.5%					±1.5%	±1.5%	±1.5%	±1.5%	±1.5%		
Linearity with Energy for >10% of full scale (c)	±1.5%					±2%						
Damage Threshold J/cm <sup>2 (b)</sup>												
<100ns	4					4						
1µs	8					5						
300µs	30					20						
2ms	50					60						
Maximum Average Power W	200					40						
Maximum Average Power Density at Maximum Power W/cm <sup>2</sup>	120 (e)					200 (e)						
Uniformity over surface	±2% ov	er central 4	10mm			±2% ove	er central 6	0mm				
Cooling	fan (see	page 103	3 for deta	ails)		conduct	ion					
Weight kg	1.2					0.5						
Compliance	CE, Chii	na RoHS				CE, Chir	na RoHS					
Version												
Part Number	7Z029	50				7Z0295	54					
Notes: (a) Calibration accuracy at various wavelengths as specified		l wavelengtl										
here. At other wavelengths, there may be an additional error up to the				nm and 2940								
value given.						d above: ±29						
Notes: (b)		or wavelengt				. For waveler . For beam si						

Notes: (c) With the "user threshold" setting set to minimum. For other settings, the spec is for >10% of full scale or greater than twice the "user threshold", whichever is greater. The user threshold is not available with LaserStar, Nova/Orion, Pulsar, USBI and Quasar. For these meters, the threshold is set to minimum and the linearity spec is >10% of full scale. The PE-C series will only operate with Nova or Orion meters with an additional adapter Ophir P/N 7Z08272 (see page 103). The adapter can introduce up to 1% additional measurement error. The user threshold feature allows adjustment of the internal threshold up to 25% of full scale if desired to avoid false triggering in noisy environments.


For further information, see the FAQs on our Website.


Notes: (d) With the LaserStar, Pulsar, USBI, Quasar and Nova/Orion with adapter only 2 of the pulse width settings are available, the 1ms and 10ms settings.


Notes: (e) For maximum power. For lower powers the damage threshold is correspondingly higher.


Notes: (f) For powers below 50W it is recommended to work with the fan off. If working with the fan on, the threshold must be set to 6% and the lowest measurable energies will be as follows:

Max Pulse Width Setting	1ms	2ms	5ms	10ms	20ms
Lowest Measurable Energy mJ	4mJ	4mJ	4mJ	4mJ	4mJ











FPE80BF-DIF-C

101

# **1.2.4** Energy Sensors Accessories

### 1.2.4.1 Accessories for Pyroelectric Sensors

Fiberoptic Adapter for Pyroelectric Sensors



Oscilloscope Adapter for Pyroelectric Sensors




Heat Sink for PE-C Series Sensors




**Beam Splitter Assembly** 










Beam splitter installed - reflected beam on sensor

Beam Splitter removed - direct beam on sensor


### **Beam Splitter Specifications**

Material	UV grade fused silica							
Spectral range	0.19 - 2.2μm							
Aperture	Ø60mm							
Damage threshold for pulses	< 10ns PW	>300µs PW						
	5J/cm <sup>2</sup> >200J/cm <sup>2</sup>							
Fraction split off See graph								



Accessory	Description	Part number	er				
Heat Sink	Heat sink that screws onto rear of PE25 and PE50 series sensors and allows working at over 50% higher average powers.	7Z08267					
Scope Adapter	Plugs in between the PE sensor and power meter. 7Z11012 Provides BNC output to scope to see every pulse up to the maximum frequency of the sensor.						
Fiber Adapters	To mount fibers to sensors you need an adapter bracket ar adapter bracket selected.	nd fiber adapt	er. All fiber a	dapters are o	compatible with the		
Fiber Adapter Brackets	Mounting brackets to allow mounting fiber adapters to py	roelectric sen	sors.				
PE Sensor Family Type		Bracket P/N		Distance from fiber to detecto			
PD10-C / PD10-IR-C / PD10-pJ-C / PD10-IR-pJ-C		7Z08275		10mm			
PE50-C / PE50BF-C PE9-C / PE9-ES-C / PE10-C / PE10BF-C / PE25-C / PE25BF-C		7Z08270 7Z08269		15mm 10mm			
Fiber Adapters	Fiber adapters for mounting to above brackets	SC type	ST type	FC type	SMA type		
For all PE sensors above		7Z08227	7Z08226	7Z08229	1G01236A		
Beam Splitter Assembly	Beam Splitter Assembly to measure pulsed laser sources too energetic for direct measurement. The reading with the Beam Splitter can be calibrated by setting the laser to a lower energy that will not damage the sensor and then taking a measurement with the beam splitter and without and taking the ratio.	7Z17001					

# 1.2.4.1 Accessories for Pyroelectric Sensors - Continued



Accessory	Description	Part number			
Damage Threshold Test Plates	Test plates with same absorber coating as the sensor. For testing that laser beam is not above damage threshold	Metallic type	BF type	THz type	BB type
		7E06031A	7E06031D	7E06031F	7E06031C
Nova PE-C Adapter	The adapter plugs between the Nova D15 socket and the smart plug of the PE-C sensor to allow the Nova to operate with PE-C series sensors. See PE-C spec sheet for details.	7Z08272			
PE-C to PE Size Adapter	The newer PE-C series sensors have a Ø62mm diameter. The older PE series sensors have a Ø85mm diameter. This adapter allows using the PE-C type sensors in jigs and setups that were originally designed for PE sensors.	7Z08273			
N Polarity Power Supply/Charger AC/DC 12V 2A N-2.1x5.5	For FPE80BF-DIF-C sensor (1 unit supplied with the sensor)	7E05029			

